


CERVICAL FIXATION SYSTEM
CERVICAL PEDICULAR SCREWS

PROXiphias-C™



Posterior Cervical Pedicular Screws



SIZE	REF.CODE
Ø 3,5 x 10 mm	101.01 3510
Ø 3,5 x 12 mm	101.01 3512
Ø 3,5 x 14 mm	101.01 3514
Ø 3,5 x 16 mm	101.01 3516
Ø 3,5 x 18 mm	101.01 3518
Ø 3,5 x 20 mm	101.01 3520
Ø 3,5 x 22 mm	101.01 3522
Ø 3,5 x 24 mm	101.01 3524
Ø 3,5 x 26 mm	101.01 3526
Ø 3,5 x 28 mm	101.01 3528
Ø 3,5 x 30 mm	101.01 3530
Ø 3,5 x 32 mm	101.01 3532
Ø 3,5 x 34 mm	101.01 3534



SIZE	REF.CODE
Ø 4,0 x 10 mm	101.01 4010
Ø 4,0 x 12 mm	101.01 4012
Ø 4,0 x 14 mm	101.01 4014
Ø 4,0 x 16 mm	101.01 4016
Ø 4,0 x 18 mm	101.01 4018
Ø 4,0 x 20 mm	101.01 4020
Ø 4,0 x 22 mm	101.01 4022
Ø 4,0 x 24 mm	101.01 4024
Ø 4,0 x 26 mm	101.01 4026
Ø 4,0 x 28 mm	101.01 4028
Ø 4,0 x 30 mm	101.01 4030
Ø 4,0 x 32 mm	101.01 4032
Ø 4,0 x 34 mm	101.01 4034

Occipital Plate



SIZE	REF.CODE
Small	121.02 0110
Medium	121.02 0115
Large	121.02 0120



Ø3,0 Screws	REF.CODE	Ø3,5 Screws	REF.CODE
Ø3,5x8 mm	121.02 6052	Ø4,0x8 mm	121.02 5036
Ø3,5x10 mm	121.02 5020	Ø4,0x10 mm	121.02 5040
Ø3,5x12 mm	121.02 5024	Ø4,0x12 mm	121.02 5044
Ø3,5x14 mm	121.02 5028	Ø4,0x14 mm	121.02 5048
Ø3,5x16 mm	121.02 5032	Ø4,0x16 mm	121.02 5052



Prodorth Posterior Cervical Fixation System offers multiple screws, hooks, rods with hexagonal tips, transition rods, occipital plates and different type of connectors for a better match to patient anatomy.

- Friction head screw feature
- Offers wide range of angulation
- Five points of fixation on the occipital plate
- Colour-coded implants for easy identification and user-friendly instrumentation to facilitate the procedure
- Self-tapping Screws for a proper and safe placement and are available in 3.5 mm & 4.0 mm diameters

Cervical Connectors

Cervical Multiaxial Connectors

PROLOBSTER™



SIZE	REF.CODE
20 mm - 30 mm	105.11 2030
30 mm - 40 mm	105.11 3040
40 mm - 50 mm	105.11 4050

Cervical Linear Connectors

PROLOBSTER-L™



SIZE	REF.CODE
30 mm	105.12 0030
40 mm	105.12 0040
50 mm	105.12 0050
60 mm	105.12 0060

PROWALRUS™



Laminar Narrow Blade Laminar Wide Blade Pedicular

SIZE	REF.CODE
Narrow Blade - Small	151.01 0017N-CS
Narrow Blade - Medium	151.01 0017N-CM
Narrow Blade - Large	151.01 0017N-CL
Wide Blade - Small	151.01 0017W-CS
Wide Blade - Medium	151.01 0017W-CM
Wide Blade - Large	151.01 0017W-CL
Pedicular - Small	151.01 0018-CS
Pedicular - Medium	151.01 0018-CM
Pedicular - Large	151.01 0018-CL

PROLOBSTER-B™



Bar Type Multiaxial Connectors

SIZE	REF.CODE
20 mm - 30 mm	105.01 2030-C
30 mm - 40 mm	105.01 3040-C

Cervical Rods


Titanium Rod

SIZE	REF.CODE
Ø 3,5 x 30 mm	104.01 3530
Ø 3,5 x 40 mm	104.01 3540
Ø 3,5 x 50 mm	104.01 3550
Ø 3,5 x 60 mm	104.01 3560
Ø 3,5 x 70 mm	104.01 3570
Ø 3,5 x 80 mm	104.01 3580
Ø 3,5 x 90 mm	104.01 3590
Ø 3,5 x 100 mm	104.01 35100
Ø 3,5 x 110 mm	104.01 35110
Ø 3,5 x 120 mm	104.01 35120
Ø 3,5 x 130 mm	104.01 35130
Ø 3,5 x 140 mm	104.01 35140
Ø 3,5 x 150 mm	104.01 35150
Ø 3,5 x 160 mm	104.01 35160
Ø 3,5 x 170 mm	104.01 35170
Ø 3,5 x 180 mm	104.01 35180
Ø 3,5 x 190 mm	104.01 35190
Ø 3,5 x 200 mm	104.01 35200




Transition Rod

SIZE	REF.CODE
Ø3,5x100 - Ø5,5x250 mm	104. 01 56350




Friction head screw feature



Offers wide range of angulation

Colour-coded implants for easy identification and user-friendly instrumentation to facilitate the procedure



 **PRODORTH**

